





# LENTIS<sup>®</sup>

## M<sub>plus</sub>

The premium IOL family for optimal treatment of cataract, presbyopia and astigmatism



PREMIUM WITHOUT COMPROMISE



DON'T  
TO BE  
FREE

# LENTIS® MplusFamily

The premium multifocal IOL family with the varifocal effect - as individual as your patients




Vision without reading glasses

Illustrative image to simulate expected and potential outcome

## Presbyopia

For many patients, the age-related deterioration of near and intermediate vision goes hand in hand with a sense of losing their freedom and independence. Multifocal intraocular lenses are an alternative to reading and varifocal spectacles and offer a permanent solution to this problem.

## Individual visual requirements need individual solutions!




Near

Intermediate

Distance

- Life style preference for intermediate ranges (computer work, driving, etc.)




Near

Intermediate

Distance

- Allrounder for active daily life style requirements



Near

Intermediate

Distance

- "X" -tra power for reading



- Large distance optic zone
- Smooth transition from distance to near
- Very high light transmission for excellent vision and contrast



Patient perception with the LENTIS® Mplus\*

## Clinical results



### Optimised vision with LENTIS® Mplus MF20


in 198 eyes

<https://www.springermedizin.de/factors-for-good-near-and-distance-visual-outcomes-of-multifocal/18031582>

**Jung Wan Kim (MD)**



- Results show excellent distance visual acuity and good near visual acuity at 40cm
- Excellent defocus performance
- Good IOL performance for younger patients with a low level of residual refraction



### Reading performance of LENTIS® Mplus MF30 and MF30T


in 34 eyes

<https://www.healio.com/ophthalmology/journals/jrs/>

**Prof. Gerd Auffarth (MD)**



- Good near and intermediate vision as well as good reading speed
- Individually preferred reading distance: 39cm
- Individually preferred intermediate distance: 62cm



### Intermediate and near vision with LENTIS® MplusX


in 50 eyes

[https://www.ajo.com/article/S0002-9394\(18\)30171-5/fulltext](https://www.ajo.com/article/S0002-9394(18)30171-5/fulltext)

**Prof. Thomas Kohnen (MD)**




- Continuously good visual acuity at all distances, especially for intermediate and near
- Good reading acuity
- Good contrast sensitivity
- High spectacle independence




# LENTIS<sup>®</sup> Mplus<sup>toric</sup> Family

## Multifocal-toric premium IOLs for correcting presbyopia and astigmatism



### Astigmatism

Impaired vision at all distances, due to a not evenly round cornea.



## Clinical results



### Results of LENTIS<sup>®</sup> Mplus | Mplus<sup>toric</sup> MF30T

in 9366 | 89 eyes

<https://www.sciencedirect.com/science/article/abs/pii/S0886335013001120>

**Jan Venter (MD)**



- Effective correction of corneal astigmatism
- Stable positioning in the capsular bag, no significant IOL rotation
- Excellent visual acuity and high patient satisfaction



### Refractive stability of the LENTIS<sup>®</sup> Mplus<sup>toric</sup> MF30T

in 70 eyes

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4886883/>

**Patrick J. Chiam (MD)**



- High rotational stability (<2.8°)
- Effective correction of the corneal cylinder
- Pre-op astigmatism average cyl 1.41 D
- Post-op residual astigmatism average cyl 0.18 D

One-piece multifocal and multifocal toric posterior chamber lenses for optimized depth of field with aspherical surface for capsular bag fixation

Parameters	LENTIS® Mplus <sup>×</sup>   Mplus	LENTIS® Mplus <sup>×toric</sup>   Mplus <sup>toric</sup>
<b>Type</b>	Foldable one-piece multifocal acrylic IOL for capsular bag fixation	Foldable one-piece multifocal <b>toric</b> acrylic IOL for capsular bag fixation
<b>Optic Size</b>	6.0 mm	6.0 mm
<b>Overall Length</b>	11.0 mm	11.0 mm
<b>Haptic Angulation</b>	0°	0°
<b>Optic Design</b>	- Dioptres: Convex-concave + Dioptres: Biconvex Aspherical surface - posterior, sector-shaped nearvision segment - anterior	Biconvex Aspherical and toric surface - posterior, sector-shaped nearvision segment- anterior Additionally available with <b>violet light filter</b>
<b>IOL Design</b>	Optics and haptics with square edges, posterior 360° continuous barrier effect	Plate haptic Optics and haptics with square edges
<b>Material</b>	HydroSmart® - a copolymer, consisting of acrylates with hydrophobic properties, UV absorbing	
<b>Available Diopters</b>	SE -10.0 D to -1.0 D (1.0 D) SE ±0.0 D to +36.0 D (0.5 D)	sph. ±0.0 D to +36.0 D (0.01 D) cyl. +0.5 D to +10.0 D (0.01 D) (sph. + cyl. < 40.0 D) axis (1°-scaling)
<b>Refractive Index</b>	1.46	1.46
<b>A constant (nominal)</b>	118.0	118.0
<b>Sterilisation</b>	Steam sterilisation	Steam sterilisation
<b>Storage</b>	Supplied in sterile water	Supplied in sterile water
<b>Recommended Injector-Sets</b> [disposable]	Check compatibility of IOL with injector matrix provided at <a href="https://lentis-eifu.com">https://lentis-eifu.com</a>	

Source: IOLcon.org

Please note that neither Teleon nor IOLcon can be held responsible for correctly specifying the optimized A constants for the Zeiss IOLMaster. The specified constants are therefore to be seen as a guide value and starting point for calculating the IOL refractive power

## Advantages of the LENTIS® Mplusfamily compared to the diffractive MIOL restoration:

- Continuous vision over all distances due to the varifocal effect
- Very high use of light
- Unique brilliant optics without rings
- Very little photic phenomena

MANUFACTURER:

Revision: QF2322v4



Teleon Surgical B.V. | Van Rensselaerweg 4 b | NL - 6956AV Spankeren

✉ [marketing@teleon-surgical.com](mailto:marketing@teleon-surgical.com)