




OPHTHALMIC SUTURE


SPECIFICATIONS

POLYPROPYLENE, BLUE MONOFILAMENT


Advanced Micro Spatula Needle

MODEL-DOUBLE ARMED			CURVATURE	WIRE	LENGTH
9-0 straight	P90062 (20cm)*		ST000-6 straight 	0°	0.15mm 16.0mm
9-0, straight and 1/2 circle	P9006062 (30cm)*		ST000-6 straight and ST160-6 1/2 circle 	160°	0.15mm 16.0mm and 6.2m
9-0 straight and 1/4 circle	P9009062*		ST000-6 straight and ST90-6 1/4 circle 	0°	0.16mm 16.0mm

Reverse Cutting

MODEL-SINGLE ARMED			CURVATURE	WIRE	LENGTH
5-0, 3/8 circle		P540151RC (70cm)*	RC-130-15 3/8 circle 	140°	0.38mm 13.1mm

* FDA 510(k) only



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA TEL: 608-270-3880 FAX: 608-270-3882 www.visionarymedicalsupplies.com