









OPHTHALMIC SUTURE


SPECIFICATIONS

NYLON, BLACK MONOFILAMENT


Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
10-0, 1/2 circle	N106062 (30cm)	N106061 (15cm)	ST160-6 1/2 circle (TG160-6)	 160°	0.15mm	5.60mm
10-0, 3/8 circle	N104062 (30cm)	N104061 (15cm)	ST140-6 3/8 circle; (TG140-6)	 140°	0.15mm	6.00mm
10-0, 3/8 circle	N104082 (40cm)	N104081 (40cm)	ST140-8 3/8 circle; (TG140-8)	 140°	0.20mm	6.00mm
9-0, 1/2 circle	N96062 (30cm)	N96061 (15cm)	ST160-6 1/2 circle (TG160-6)	 160°	0.15mm	5.60mm
9-0, 3/8 circle	N94062 (30cm)	N94061 (15cm)	ST140-6 3/8 circle (TG140-6)	 140°	0.15mm	6.00mm
9-0, 3/8 circle	N94082 (40cm)	N94081 (40cm)	ST140-8 3/8 circle (TG140-8)	 140°	0.20mm	6.00mm
8-0, 3/8 circle	N84082 (30cm)*		ST140-8 3/8 circle (TG140-8)	 140°	0.20mm	6.00mm
6-0, 3/8 circle	N640112 (76cm)*		ST140-11 3/8 circle	 140°	0.28mm	8.00mm

Taper Point

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
7-0, 3/8 circle	N74092 (30cm)*		TP140-9 3/8 circle (BV130-5)	 140°	0.23mm	6.00mm

* FDA 510(k) only



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE