



OPHTHALMIC SUTURE


SPECIFICATIONS

SILK, BLACK BRAIDED


Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	S84082 (30cm)	S84081 (15cm)	ST140-8 3/8 circle (TG140-8)		140°	0.20mm 6.00mm
6-0, 3/8 circle	S640112 (76cm)*		ST140-11 3/8 circle		140°	0.28mm 8.00mm

Taper Point


	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
7-0, 3/8 circle	S74092 (45cm)*		TP140-9 3/8 circle (BV130-5)		140°	0.23mm 6.00mm

Cutting


	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
4-0, 3/8 circle		S440171 (76cm)*	CT-135-17 3/8 circle		140°	0.23mm 6.00mm

SILK, VIRGIN

Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	VS84082 (40cm)	VS84081 (40cm)	ST140-8 3/8 circle (TG140-8)		140°	0.20mm 6.00mm

* FDA 510(k) only



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE