



OPHTHALMED™

LASER PROBES

FOR ENDO-OCULAR
PHOTOCOAGULATION
TREATMENTS

WWW.OPHTHALMED.COM

(877) 4-25 GAUGE (4284)


In business since 2001, Ophthalmed LLC specializes in designing and manufacturing higher quality endo photocoagulation laser probes.

The range of our products is one of the most extensive on the market, backed up by several patent pending applications, and private labeled to major companies in US and worldwide, including well known laser machine manufacturers and other high profile distributors.

We are committed to serving the needs of the ophthalmic industry, and to providing the latest technologies to the surgeons, including custom designs.

Our current product range includes laser and illuminated laser probes in various tip sizes and configurations, including extended (patent pending). We connect to most laser and illumination machines.

Our probes are manufactured in accordance with a quality system certified to ISO 9001 and ISO 13485 CMDCAS, and cleared/approved by various regulatory agencies, such as US FDA, Health Canada and the European Community.



ADDITIONAL VARIATIONS ARE AVAILABLE.

ALCON / COHERENT / ZEISS / ELLEX CONNECTOR

A217000	Straight	
A237000	Angled	20g
A257000	Curved	
A827000	Straight	23g
A877000	Angled	
A327000	Straight	
A377000	Angled	25g
A367000	Curved	
AA27000	Straight	27g
AA67000	Curved	



NIDEK CONNECTOR - THREADED

A212000	Straight	
A232000	Angled	20g
A252000	Curved	
A822000	Straight	23g
A872000	Angled	
A322000	Straight	
A372000	Angled	25g
A362000	Curved	
AA22000	Straight	27g
AA62000	Curved	



IRIDEX CONNECTOR

A218000	Straight	
A238000	Angled	20g
A258000	Curved	
A828000	Straight	23g
A878000	Angled	
A328000	Straight	
A378000	Angled	25g
A368000	Curved	
AA28000	Straight	27g
AA68000	Curved	



NIDEK CONNECTOR - PUSH-IN

A215000	Straight	
A235000	Angled	20g
A255000	Curved	
A825000	Straight	23g
A875000	Angled	
A325000	Straight	
A375000	Angled	25g
A365000	Curved	
AA25000	Straight	27g
AA65000	Curved	



QUANTEL VITRA / SUPRA CONNECTOR (DIRECT)

A219000	Straight	
A239000	Angled	20g
A259000	Curved	
A829000	Straight	23g
A879000	Angled	
A329000	Straight	
A379000	Angled	25g
A369000	Curved	
AA29000	Straight	27g
AA69000	Curved	




HGM CONNECTOR - NO BAR CODE

A214000	Straight	
A234000	Angled	20g
A254000	Curved	
A824000	Straight	23g
A874000	Angled	
A324000	Straight	
A374000	Angled	25g
A364000	Curved	
AA24000	Straight	27g
AA64000	Curved	



STANDARD NON-EXTENDED TIPS



EXTENDED TIPS - TYPICAL



PROTECTED BY ONE OR MORE PATENT PENDING APPLICATIONS.
ADDITIONAL VARIATIONS ARE AVAILABLE.

ALCON / COHERENT / ZEISS / ELLEX LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A527200	Straight	
A547200	Angled	
A567200	Curved	20g
A587200	Straight, Extended	
A597200	Angled, Extended	
A5A7200	Curved, Extended	
B527203	Straight	
B577203	Angled	23g
B587203	Straight, Extended	
B5B7203	Angled, Extended	
B527204	Straight	
B567204	Curved	25g
B587204	Straight, Extended	
B5A7204	Curved, Extended	
B527205	Straight	27g
B567205	Curved	



PROTECTED BY ONE OR MORE
PATENT PENDING APPLICATIONS.

ADDITIONAL VARIATIONS ARE AVAILABLE.

19 GAUGE ILLUMINATION CONNECTOR

A527100	Straight	
A547100	Angled	
A567100	Curved	20g
A587100	Straight, Extended	
A597100	Angled, Extended	
A5A7100	Curved, Extended	
B527103	Straight	
B577103	Angled	23g
B587103	Straight, Extended	
B5B7103	Angled, Extended	
B527104	Straight	
B567104	Curved	25g
B587104	Straight, Extended	
B5A7104	Curved, Extended	
B527105	Straight	27g
B567105	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A527300	Straight	
A547300	Angled	
A567300	Curved	20g
A587300	Straight, Extended	
A597300	Angled, Extended	
A5A7300	Curved, Extended	
B527303	Straight	
B577303	Angled	23g
B587303	Straight, Extended	
B5B7303	Angled, Extended	
B527304	Straight	
B567304	Curved	25g
B587304	Straight, Extended	
B5A7304	Curved, Extended	
B527305	Straight	27g
B567305	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A527500	Straight	
A547500	Angled	
A567500	Curved	20g
A587500	Straight, Extended	
A597500	Angled, Extended	
A5A7500	Curved, Extended	
B527503	Straight	
B577503	Angled	23g
B587503	Straight, Extended	
B5B7503	Angled, Extended	
B527504	Straight	
B567504	Curved	25g
B587504	Straight, Extended	
B5A7504	Curved, Extended	
B527505	Straight	27g
B567505	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A527400	Straight	
A547400	Angled	
A567400	Curved	20g
A587400	Straight, Extended	
A597400	Angled, Extended	
A5A7400	Curved, Extended	
B527403	Straight	
B577403	Angled	23g
B587403	Straight, Extended	
B5B7403	Angled, Extended	
B527404	Straight	
B567404	Curved	25g
B587404	Straight, Extended	
B5A7404	Curved, Extended	
B527405	Straight	27g
B567405	Curved	



IRIDEX LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A528200	Straight	
A548200	Angled	
A568200	Curved	20g
A588200	Straight, Extended	
A598200	Angled, Extended	
A5A8200	Curved, Extended	
<hr/>		
B528203	Straight	
B578203	Angled	23g
B588203	Straight, Extended	
B5B8203	Angled, Extended	
<hr/>		
B528204	Straight	
B568204	Curved	25g
B588204	Straight, Extended	
B5A8204	Curved, Extended	
<hr/>		
B528205	Straight	27g
B568205	Curved	



PROTECTED BY ONE OR MORE PATENT PENDING APPLICATIONS.

ADDITIONAL VARIATIONS ARE AVAILABLE.

19 GAUGE ILLUMINATION CONNECTOR

A528100	Straight	
A548100	Angled	
A568100	Curved	20g
A588100	Straight, Extended	
A598100	Angled, Extended	
A5A8100	Curved, Extended	
<hr/>		
B528103	Straight	
B578103	Angled	23g
B588103	Straight, Extended	
B5B8103	Angled, Extended	
<hr/>		
B528104	Straight	
B568104	Curved	25g
B588104	Straight, Extended	
B5A8104	Curved, Extended	
<hr/>		
B528105	Straight	27g
B568105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A528500	Straight	
A548500	Angled	
A568500	Curved	20g
A588500	Straight, Extended	
A598500	Angled, Extended	
A5A8500	Curved, Extended	
<hr/>		
B528503	Straight	
B578503	Angled	23g
B588503	Straight, Extended	
B5B8503	Angled, Extended	
<hr/>		
B528504	Straight	
B568504	Curved	25g
B588504	Straight, Extended	
B5A8504	Curved, Extended	
<hr/>		
B528505	Straight	27g
B568505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A528300	Straight	
A548300	Angled	
A568300	Curved	20g
A588300	Straight, Extended	
A598300	Angled, Extended	
A5A8300	Curved, Extended	
<hr/>		
B528303	Straight	
B578303	Angled	23g
B588303	Straight, Extended	
B5B8303	Angled, Extended	
<hr/>		
B528304	Straight	
B568304	Curved	25g
B588304	Straight, Extended	
B5A8304	Curved, Extended	
<hr/>		
B528305	Straight	27g
B568305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A528400	Straight	
A548400	Angled	
A568400	Curved	20g
A588400	Straight, Extended	
A598400	Angled, Extended	
A5A8400	Curved, Extended	
<hr/>		
B528403	Straight	
B578403	Angled	23g
B588403	Straight, Extended	
B5B8403	Angled, Extended	
<hr/>		
B528404	Straight	
B568404	Curved	25g
B588404	Straight, Extended	
B5A8404	Curved, Extended	
<hr/>		
B528405	Straight	27g
B568405	Curved	



NIDEK LASER CONNECTOR - THREADED

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A522200	Straight	
A542200	Angled	
A562200	Curved	20g
A582200	Straight, Extended	
A592200	Angled, Extended	
A5A2200	Curved, Extended	

B522203	Straight	
B572203	Angled	23g
B582203	Straight, Extended	
B5B2203	Angled, Extended	

B522204	Straight	
B562204	Curved	25g
B582204	Straight, Extended	
B5A2204	Curved, Extended	

B522205	Straight	27g
B562205	Curved	



PROTECTED BY ONE OR MORE
PATENT PENDING APPLICATIONS.

ADDITIONAL VARIATIONS ARE AVAILABLE.

19 GAUGE ILLUMINATION CONNECTOR

A522100	Straight	
A542100	Angled	
A562100	Curved	20g
A582100	Straight, Extended	
A592100	Angled, Extended	
A5A2100	Curved, Extended	

B522103	Straight	
B572103	Angled	23g
B582103	Straight, Extended	
B5B2103	Angled, Extended	

B522104	Straight	
B562104	Curved	25g
B582104	Straight, Extended	
B5A2104	Curved, Extended	

B522105	Straight	27g
B562105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A522500	Straight	
A542500	Angled	
A562500	Curved	20g
A582500	Straight, Extended	
A592500	Angled, Extended	
A5A2500	Curved, Extended	

B522503	Straight	
B572503	Angled	23g
B582503	Straight, Extended	
B5B2503	Angled, Extended	

B522504	Straight	
B562504	Curved	25g
B582504	Straight, Extended	
B5A2504	Curved, Extended	

B522505	Straight	27g
B562505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A522300	Straight	
A542300	Angled	
A562300	Curved	20g
A582300	Straight, Extended	
A592300	Angled, Extended	
A5A2300	Curved, Extended	

B522303	Straight	
B572303	Angled	23g
B582303	Straight, Extended	
B5B2303	Angled, Extended	

B522304	Straight	
B562304	Curved	25g
B582304	Straight, Extended	
B5A2304	Curved, Extended	

B522305	Straight	27g
B562305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A522400	Straight	
A542400	Angled	
A562400	Curved	20g
A582400	Straight, Extended	
A592400	Angled, Extended	
A5A2400	Curved, Extended	

B522403	Straight	
B572403	Angled	23g
B582403	Straight, Extended	
B5B2403	Angled, Extended	

B522404	Straight	
B562404	Curved	25g
B582404	Straight, Extended	
B5A2404	Curved, Extended	

B522405	Straight	27g
B562405	Curved	



QUANTEL VITRA / SUPRA LASER CONNECTOR (DIRECT)

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A529200	Straight	
A549200	Angled	
A569200	Curved	20g
A589200	Straight, Extended	
A599200	Angled, Extended	
A5A9200	Curved, Extended	
B529203	Straight	
B579203	Angled	23g
B589203	Straight, Extended	
B5B9203	Angled, Extended	
B529204	Straight	
B569204	Curved	25g
B589204	Straight, Extended	
B5A9204	Curved, Extended	
B529205	Straight	27g
B569205	Curved	



PROTECTED BY ONE OR MORE PATENT PENDING APPLICATIONS.

ADDITIONAL VARIATIONS ARE AVAILABLE.

19 GAUGE ILLUMINATION CONNECTOR

A529100	Straight	
A549100	Angled	
A569100	Curved	20g
A589100	Straight, Extended	
A599100	Angled, Extended	
A5A9100	Curved, Extended	
B529103	Straight	
B579103	Angled	23g
B589103	Straight, Extended	
B5B9103	Angled, Extended	
B529104	Straight	
B569104	Curved	25g
B589104	Straight, Extended	
B5A9104	Curved, Extended	
B529105	Straight	27g
B569105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A529500	Straight	
A549500	Angled	
A569500	Curved	20g
A589500	Straight, Extended	
A599500	Angled, Extended	
A5A9500	Curved, Extended	
B529503	Straight	
B579503	Angled	23g
B589503	Straight, Extended	
B5B9503	Angled, Extended	
B529504	Straight	
B569504	Curved	25g
B589504	Straight, Extended	
B5A9504	Curved, Extended	
B529505	Straight	27g
B569505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A529300	Straight	
A549300	Angled	
A569300	Curved	20g
A589300	Straight, Extended	
A599300	Angled, Extended	
A5A9300	Curved, Extended	
B529303	Straight	
B579303	Angled	23g
B589303	Straight, Extended	
B5B9303	Angled, Extended	
B529304	Straight	
B569304	Curved	25g
B589304	Straight, Extended	
B5A9304	Curved, Extended	
B529305	Straight	27g
B569305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR


A529400	Straight	
A549400	Angled	
A569400	Curved	20g
A589400	Straight, Extended	
A599400	Angled, Extended	
A5A9400	Curved, Extended	
B529403	Straight	
B579403	Angled	23g
B589403	Straight, Extended	
B5B9403	Angled, Extended	
B529404	Straight	
B569404	Curved	25g
B589404	Straight, Extended	
B5A9404	Curved, Extended	
B529405	Straight	27g
B569405	Curved	



NIDEK LASER CONNECTOR - PUSH-IN

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A525200	Straight	
A545200	Angled	
A565200	Curved	20g
A585200	Straight, Extended	
A595200	Angled, Extended	
A5A5200	Curved, Extended	
<hr/>		
B525203	Straight	
B575203	Angled	23g
B585203	Straight, Extended	
B5B5203	Angled, Extended	
<hr/>		
B525204	Straight	
B565204	Curved	25g
B585204	Straight, Extended	
B5A5204	Curved, Extended	



PROTECTED BY ONE OR MORE
PATENT PENDING APPLICATIONS.

ADDITIONAL VARIATIONS ARE AVAILABLE.

19 GAUGE ILLUMINATION CONNECTOR

A525100	Straight	
A545100	Angled	
A565100	Curved	
A585100	Straight, Extended	20g
A595100	Angled, Extended	
A5A5100	Curved, Extended	
<hr/>		
B525103	Straight	
B575103	Angled	23g
B585103	Straight, Extended	
B5B5103	Angled, Extended	
<hr/>		
B525104	Straight	
B565104	Curved	25g
B585104	Straight, Extended	
B5A5104	Curved, Extended	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A525500	Straight	
A545500	Angled	
A565500	Curved	
A585500	Straight, Extended	20g
A595500	Angled, Extended	
A5A5500	Curved, Extended	
<hr/>		
B525503	Straight	
B575503	Angled	23g
B585503	Straight, Extended	
B5B5503	Angled, Extended	
<hr/>		
B525504	Straight	
B565504	Curved	25g
B585504	Straight, Extended	
B5A5504	Curved, Extended	



ALCON OCUTOME ILLUMINATION CONNECTOR

A525300	Straight	
A545300	Angled	
A565300	Curved	
A585300	Straight, Extended	20g
A595300	Angled, Extended	
A5A5300	Curved, Extended	
<hr/>		
B525303	Straight	
B575303	Angled	23g
B585303	Straight, Extended	
B5B5303	Angled, Extended	
<hr/>		
B525304	Straight	
B565304	Curved	25g
B585304	Straight, Extended	
B5A5304	Curved, Extended	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A525400	Straight	
A545400	Angled	
A565400	Curved	
A585400	Straight, Extended	20g
A595400	Angled, Extended	
A5A5400	Curved, Extended	
<hr/>		
B525403	Straight	
B575403	Angled	23g
B585403	Straight, Extended	
B5B5403	Angled, Extended	
<hr/>		
B525404	Straight	
B565404	Curved	25g
B585404	Straight, Extended	
B5A5404	Curved, Extended	



HGM LASER CONNECTOR - NO BAR CODE

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A524200	Straight	
A544200	Angled	
A564200	Curved	20g
A584200	Straight, Extended	
A594200	Angled, Extended	
A5A4200	Curved, Extended	
B524203	Straight	
B574203	Angled	23g
B584203	Straight, Extended	
B5B4203	Angled, Extended	
B524204	Straight	
B564204	Curved	25g
B584204	Straight, Extended	
B5A4204	Curved, Extended	



PROTECTED BY ONE OR MORE PATENT PENDING APPLICATIONS.

ADDITIONAL VARIATIONS ARE AVAILABLE.

19 GAUGE ILLUMINATION CONNECTOR

A524100	Straight	
A544100	Angled	
A564100	Curved	20g
A584100	Straight, Extended	
A594100	Angled, Extended	
A5A4100	Curved, Extended	
B524103	Straight	
B574103	Angled	23g
B584103	Straight, Extended	
B5B4103	Angled, Extended	
B524104	Straight	
B564104	Curved	25g
B584104	Straight, Extended	
B5A4104	Curved, Extended	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A524500	Straight	
A544500	Angled	
A564500	Curved	20g
A584500	Straight, Extended	
A594500	Angled, Extended	
A5A4500	Curved, Extended	
B524503	Straight	
B574503	Angled	23g
B584503	Straight, Extended	
B5B4503	Angled, Extended	
B524504	Straight	
B564504	Curved	25g
B584504	Straight, Extended	
B5A4504	Curved, Extended	



ALCON OCUTOME ILLUMINATION CONNECTOR

A524300	Straight	
A544300	Angled	
A564300	Curved	20g
A584300	Straight, Extended	
A594300	Angled, Extended	
A5A4300	Curved, Extended	
B524303	Straight	
B574303	Angled	23g
B584303	Straight, Extended	
B5B4303	Angled, Extended	
B524304	Straight	
B564304	Curved	25g
B584304	Straight, Extended	
B5A4304	Curved, Extended	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A524400	Straight	
A544400	Angled	
A564400	Curved	20g
A584400	Straight, Extended	
A594400	Angled, Extended	
A5A4400	Curved, Extended	
B524403	Straight	
B574403	Angled	23g
B584403	Straight, Extended	
B5B4403	Angled, Extended	
B524404	Straight	
B564404	Curved	25g
B584404	Straight, Extended	
B5A4404	Curved, Extended	





Aspiration add-on functionality is available upon request.

A flexible plastic jacket covers and protects the 8ft long fiber.

Homogeneous energy distribution as well as round and consistently small laser spot allow the tip to stay few millimeters away from the retina.

Ergonomic and non-slip handpiece design.



Different connectors that suit most laser and illumination machines.

Sterilized by Ethylene Oxide, double pouched, single-use, latex-free.

Available in Dia 0.41 mm (27ga), Dia 0.51 mm (25ga), Dia 0.64 mm (23ga), and Dia 0.91 mm (20g). Length 2.4 m (8 ft). Box of 10.

Do not re-sterilize. Do not reuse.

Material fatigue decreases performance.

WARNING

Follow the recommendations of the laser and illumination source manufacturers (e.g.: calibration of the laser source as well as the need for protective eye wear for the corresponding wavelength of the laser model used).

WARNING IN CASE OF ILLUMINATION

Due to the potentially hazardous effects of phototoxicity to the retina from prolonged endoillumination, this device should be used with an illumination source with filters that eliminate UV radiation (<400nm) and whenever possible, with filters that eliminate short wavelength blue light (<420nm). Care should be taken not to focus this device on a single point of the retina for prolonged, unnecessary periods of time.

CAUTION: FEDERAL (U.S.A) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

CONTACT INFORMATION:

+1 770 777 6613

+1 877 4 25 GAUGE (4284)

+1 678 623 3765 (FAX)

info@ophthalmed.com

1050 Northfield Court, Suite 280
Roswell, GA 30076 USA

DISTRIBUTED BY:



OPHTHALMEDTM

MADE IN THE USA