

Informações Biometro Optico

	IOL	Nominal	Haigis	HofferQ	Holl.1	SRT/T	SRK II	*Barret	Holl.2	*Hill RBF	*KANE
	LENTIS Comfort LS-313 MF15	A = 118.0	a0 = 1.019 a1 = 0.309 a2 = 0.107	pACD = 5.15	sf = 1.38	A = 118.32	A = 118.5	LF = 1.53	5.150	A = 118.32	A = 118.32
	LENTIS Mplus LS-313 MF20	A=118.0	a0 = 1.142 a1 = 0.400 a2 = 0.100	pACD = 5.38	sf = 1.57	A = 118.66	A = 118.9	LF = 1.71 DF = 0	5.350	A = 118.66	A = 118.66
	LENTIS Mplus LS-313 MF30	A = 118.0	a0 = 0.950 a1 = 0.400 a2 = 0.100	pACD = 5.21	sf = 1.47	A = 118.50	A = 118.6	LF = 1.62 DF = 0	5.260	A = 118.50	A = 118.50
	FEMTIS Comfort FB-313 MF15	A = 117.8	a0 = 0.759 a1 = 0.400 a2 = 0.100	pACD = 5.03	sf = 1.27	A = 118.14	A = 118.4	LF = 1.43 DF = 0	5.050	A = 118.14	A = 118.14